



Fingerprint Development Chamber


- Detect Quality Latent Prints on Porous Surfaces with Reproducible Results



SD-34S
with optional Vent-Box



Product



34S



Schedule
Contract # GS-47QSMS24D0042



CONTENTS:

- Product Overview (p.2)
- Design Features (p.3)
- Performance & Selection (p.4)
- Filtration Technology (p.5)
- Specifications (p.6)
- Options & Accessories (p.8)

SafeDEVELOP™

Fingerprint Development Chamber

34S

PRODUCT OVERVIEW

2

INTRODUCTION

Safe DEVELOP™ Fingerprint Development Chambers are designed to accelerate the processing of latent fingerprints on porous surfaces using DFO, Ninhydrin and other development chemicals within a controlled environment for optimum effectiveness where moisture, temperature and time are critical factors.

Safe DEVELOP Fingerprint Development Chamber controls all functions from start-to-finish, permitting the investigator to initiate an unattended cycle, establish the proper development intensity and duration, and to return upon completion to collect results. The resulting prints will fluoresce with the use of various lasers and light sources.



Deep into its second generation, Air Science embraces the diversity and cultural heritage of the founders and co-workers who are continuing a tradition of excellence. Demonstrating a commitment to adaptation, inclusion, and quality output from a United States-based company with a domestic and global reach.

APPLICATIONS

Using innovative filtration technology, the Safe DEVELOP Fingerprint Development Chamber creates a safe work environment over the widest range of applications in the industry.

State and Federal Crime Laboratories \ Crime Scene Investigation \ Law Enforcement Agencies \ Medical Examiners' Programs \ Criminal Justice Education

KEY FEATURES

- The professionally designed automatic development chamber eliminates reliance on do-it-yourself systems that lack user-safety allowances and controls needed for repeatable results.
- To change between profiles requires only the press of a few buttons on an easy to use LED display.
- The chamber features rapid condition recovery following a door opening.
- The wide temperature and humidity range allows fingerprints to be processed in a matter of minutes - not days like some conventional methods.
- The steam generator produces humidity by a sealed, automatic, steam injection system, ensuring only vaporized water enters the sample chamber to eliminate the possibility of contaminating samples with water droplets.
- The chamber's generous working area and flexible shelving and rod system enable large batches of operational material to be processed quickly and easily.

SAFE DEVELOP TECHNOLOGY

DFO and Ninhydrin fuming are the most effective techniques for detecting latent prints on paper and similar porous surfaces. Safe DEVELOP performs well with DFO and Ninhydrin, however is not limited in scope and can also utilize other chemical developers, allowing illumination with various lasers and light sources to make the resulting prints fluoresce.

Development Chamber Process Table on [page 6](#).



Safedevelop SD-34S

CONTENTS:

- Product Overview (p.2)
- Design Features (p.3)
- Performance & Selection (p.4)
- Filtration Technology (p.5)
- Specifications (p.6)
- Options & Accessories (p.8)

DESIGN FEATURES

- A. Door Key:** Chamber access keys prevent unauthorized removal of evidence or accidental operator exposure to chemical fumes or high heat.
- B. Control Panel:** Front-mounted control panel with electronic On/Off switch, lights, temperature and humidity controller, cycle complete lamp and low water alarm lamp.
- C. Glass Door:** Multi-pane, heated glass door minimizes condensation and activates automatically whenever humidity is used.
- D. Low Water Light:** Low water level warning light notifies the operator when to add water.
- E. Push-Pull™ Shelving:** Perforated Push-Pull™ shelves can slide in or out with one hand.
- F. Lighting:** Dual LED light strips on the door provide a large illumination area to watch fingerprints being developed.
- G. Hanging Rods:** Removable stainless steel hanging rods with clips.
- H. Water Carboy:** Side mounted, 1 gallon (4 liter). Easy to visualize water levels and quick release coupling for simple removal and filling.
- I. Chamber:** Corrosion-resistant insulated internal chamber.
- J. Levelling:** Adjustable leveling feet.
- K. Stand:** Optional mobile cart with locking casters.
- L. Modular Filtration:** Optional Vent-Box™ filtration unit available with Multiplex™ filtration technology, a unique configuration that includes a pre-filter and main carbon filter.

ADDITIONAL FEATURES

Quality Door Construction: The large viewing area offers easy observation of critical samples along with dual vertical LED lights. The multi-pane, heated glass door minimizes condensation.

Steam Generator: The Air Science® steam generator adds heat while humidifying, providing for quick ramp-up and rapid condition recovery after door openings. Advanced humidity sensors with built-in temperature compensation provide accurate readings at all temperatures.

Preset Profiles: Safe DEVELOP cabinets are constructed of polypropylene that does not absorb liquids, is easily cleaned with household detergents and can be sprayed with a 10% bleach solution to eliminate biological contaminants.



Safedevelop SD-34S with optional Vent-Box

CONTENTS:

- Product Overview (p.2)
- Design Features (p.3)
- Performance & Selection (p.4)
- Filtration Technology (p.5)
- Specifications (p.6)
- Options & Accessories (p.8)

Each Air Science Safe DEVELOP Fingerprint Development Chamber includes features expressed through sound design and certified quality construction. Options and accessories add functional performance to meet specific applications.

PERFORMANCE

The Air Science [Multiplex Filtration System](#) offers a range of options for high performance protection.


DESIGN

Professional quality Air Science Safe DEVELOP Fingerprint Development Chambers comply with current technical and safety regulations. The cabinet frame and work surfaces, comprised of industrial components, are durable and chemically resistant.

The Air Science filter assembly is easy to access, easy to change, plus a unique filter clamping design eliminates bypass leakage outside the cabinet.

RELIABILITY

Internal systems are isolated from fumes, extending product life.



345

SELECTION

Fingerprint development chambers are available in 2 electrically specific sizes, totaling 2 standard models.



Safedevelop SD-345 with optional Vent-Box



FSA/Autocal Control Panel with On/Off switch for unit, lights, temperature and humidity controller, cycle complete lamp and low water alarm lamp. Also incorporates profile select and emergency stop switches to ensure complete unit control.

CONTENTS:

- Product Overview (p.2)
- Design Features (p.3)
- Performance & Selection (p.4)
- Filtration Technology (p.5)
- Specifications (p.6)
- Options & Accessories (p.8)

SafeDEVELOP™

Fingerprint Development Chamber

345

FILTRATION TECHNOLOGY

5



FILTRATION

At the heart of the optional Vent-Box product line is innovative filtration technology. **The Multiplex Filtration System** consists of a pre-filter, main activated carbon.

The Air Science **carbon filtration technique** is based on enhanced, activated carbon particle formulations from specially selected, naturally occurring raw material that is superior to wood or other organic sources. The carbon is treated to attain the proper porosity and aggregate surface area and to react with several ranges of aerosolized chemicals moved through the filter by an air handling blower.

SECUR.
safe disposal service

Filter disposal services are available in selected markets providing responsible destruction or recycling of used saturated filters in authorized facilities.


FILTER CONFIGURATION

The vented chemical storage cabinet can be equipped with a single activated carbon main filter activated to adsorb one or more specific vapors or family of vapors.

The carbon filter is sized to fit the specified product model number and configured to optimize airflow across 100% of the filter surface area. The self-contained assembly maximizes filter efficiency, prolongs filter life, optimizes diffusion and saturation and improves user safety.


P. Electrostatic Pre-Filter: Protects the main filters from aerosols, mists, dust and particulates.

C. Activated Carbon Main Filter: A single filter configuration.



AIRFLOW

The optional Vent-Box filtration unit maintains a constant face velocity in compliance with USA and international standards for safety and performance, creating uniform conditions within the chamber. Contaminated air is pulled through the Multiplex Filtration System; clean air is returned to the room.





CONTENTS:

- Product Overview (p.2)
- Design Features (p.3)
- Performance & Selection (p.4)
- Filtration Technology (p.5)
- Specifications (p.6)
- Options and Accessories (p.8)

PROCESS TYPE	DEVELOPMENT CHAMBER PROCESS TABLE - PROFILE SET UP CHARACTERISTIC						
	Evidence Preparation	Temp C°	RH %	Time (min)	Print Results	Photography	Unit Profile
Ninhydrin (2,2-Dihydroxyindane-1, 3-dione)	Follow proper safety precautions. Using a fume hood, apply chemical to specimen (spray, dip, brush). Completely dry before processing in chamber.	80	65	3	Purple. Repeat process as needed.	530-555 nm light source with no filter	Ninhydrin
DFO (1,8-Diazafluoren-9-one)	Follow proper safety precautions. Using a fume hood, apply chemical to specimen (spray, dip) for 5 seconds. Completely dry before processing in chamber.	100	–	20	Yellow.	495-550 nm light source with orange filter	DFO
Nickel Nitrate	Apply only after processing specimen with Ninhydrin. Follow proper safety precautions. Using a fume hood, apply chemical to specimen (spray). Completely dry before processing in chamber.	80	65	20	Ridge detail enhancement.	Green filter or 530 nm light source with no filter	Ninhydrin
5-MTN (5-Methylthioninhydrin)	Follow proper safety precautions. Using a fume hood, apply chemical to specimen (spray, dip) for 5 seconds. Completely dry before processing in chamber.	80	65	3	Strong purple, repeat if needed or try with Nickel Nitrate. Repeat process as needed.	Green filter	Ninhydrin
1,2-Indanedione	Follow proper safety precautions. Using a fume hood, apply chemical to specimen (spray, dip, wash). Completely dry before processing in chamber.	100	–	10	Pale pink. Repeat process as needed.	515 nm light source with orange filter	DFO
Zinc Chloride	Apply only after using Ninhydrin or 5-MTN to enhance prints. Follow proper safety precautions. Using a fume hood, apply chemical to specimen. Completely dry before processing in chamber.	80	65	20-40	Orange if prints previously treated with Ninhydrin. Red if prints previously treated with 5-MTN. Repeat process as needed.	Orange filter	Ninhydrin

CONTENTS:

- Product Overview (p.2)
- Design Features (p.3)
- Performance & Selection (p.4)
- Filtration Technology (p.5)
- Specifications (p.6)
- Options & Accessories (p.8)



MODEL	DIMENSIONS			WEIGHT (LBS/KG)		
	Internal Height	Work Space (W x D x H)	External (W x D x H)	Shipping (W x D x H)	Net	Ship
Safe DEVELOP						
SD-34S	19.5" / 495 mm	18.25" x 20.5" x 19.6" / 464 x 521 x 495 mm	36.25" x 25.5" x 29.25" / 921 x 648 x 743 mm	40" x 48" x 45" / 1016 x 1219 x 1143 mm	200 / 91	245 / 111
VB60 (optional)	-	-	21.25" x 12" x 13" / 540 x 305 x 330 mm	24" x 22" x 18" / 610 x 559 x 457 mm	32 / 14	34 / 15

CONTENTS:

- Product Overview (p.2)
- Design Features (p.3)
- Performance & Selection (p.4)
- Filtration Technology (p.5)
- Specifications (p.6)
- Options & Accessories (p.8)

SafeDEVELOP™

Fingerprint Development Chamber

34S

SPECIFICATIONS

8

PRODUCT SPECIFICATIONS

Construction

Airflow	<... Crossflow. ...>
Controls	<... Programmable heat and humidity controllers. Lights On/Off. ...>
Electrical	<... 120V, 60Hz or 230V, 50Hz voltages available. Specify when ordering. ...>
Monitoring	<... Audio and visual, timer cycle complete, low water level. ...>
Lighting	<... LED. ...>
Temperature Range	<... Up to 100°C (depends on profile). ...>
Relative Humidity	<... Up to 65% RH (depends on profile). ...>
Shelves	<... 2 Push-Pull, perforated on sliding rails. ...>
Hanging Rods	<... 4 stainless steel rods with 8 clips. ...>
Alarms (Audio and Visual)	<... Timer Cycle Complete, Low Water Level. ...>
Water Bottle	<... Carboy, water fill bottle with cap, 1 gallon (4 liters). Universal side mount holder mounts on side or top of unit. Includes all tubes and quick release fittings. ...>

Specifications are subject to change without notice

FILTER SPECIFICATIONS

Safe DEVELOP Model

34S with Optional Vent-Box Only

Primary Filter*	(1)
Pre-Filter*	(1)

* For specific examples refer to Multiplex filtration system summary on [page 5](#).



Through our partner company [Filtco Filters](#), Air Science is a single source supplier of all pre-filters, carbon filters and HEPA/ULPA filters used in our products.

CONTENTS:

- Product Overview (p.2)
- Design Features (p.3)
- Performance & Selection (p.4)
- Filtration Technology (p.5)
- Specifications (p.6)
- Options & Accessories (p.8)

OPTIONS & ACCESSORIES

Safe DEVELOP Model		SD-34S
Vent-Box Filtration	Modular ductless filtration system. Utilizes the Multiplex carbon filtration system with a pre-filter and main filter. Optional HEPA/ULPA filters are also available.	VB60
Heavy Duty Base Stand	Provides a lower storage half shelf; accommodates wheelchair access. Locking casters.	CART-30

CONTENTS:

- Product Overview (p.2)
- Design Features (p.3)
- Performance & Selection (p.4)
- Filtration Technology (p.5)
- Specifications (p.6)
- Warranty (p.8)

WARRANTY

This product is protected by the Air Science Legacy Limited Lifetime Warranty™.



For details visit the [Warranty section](#) of our website.

STANDARDS & COMPLIANCE

Quality Management Systems	ISO 9001 : 2015
Environment	ISO 14001: 2015 ENERGY STAR® Partner



120 6th Street \ Fort Myers, FL 33907
T. 239-489-0024 \ **Toll Free.** 800-306-0656 \ F. 800-306-0677
www.airscience.com

The information contained in this manual and the accompanying product are copyrighted and all rights are reserved by Air Science. Air Science reserves the right to make periodic minor design changes without obligation to notify any person or entity of such change.



Schedule
Contract # GS-47QSMS24D0042

